

For Research Use Only

MYH7 Recombinant antibody, PBS Only

Catalog Number: 85322-1-PBS



Basic Information

Catalog Number:

85322-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000257

GeneID (NCBI):

4625

UNIPROT ID:

P12883

Full Name:

myosin, heavy chain 7, cardiac muscle, beta

Calculated MW:

223 kDa

Observed MW:

200-220 kDa

Purification Method:

Protein A purification

CloneNo.:

242250F6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 85322-1-RR (MYH7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85322-1-PBS in a different storage buffer formulation.