

For Research Use Only

# DTX3L Recombinant antibody, PBS Only

Catalog Number: 85340-2-PBS



## Basic Information

### Catalog Number:

85340-2-PBS

### Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG2568

### GenBank Accession Number:

BC060509

### GeneID (NCBI):

151636

### UNIPROT ID:

Q8TDB6

### Full Name:

deltex 3-like (Drosophila)

### Calculated MW:

740 aa, 84 kDa

### Observed MW:

84 kDa

### Purification Method:

Protein A purification

### CloneNo.:

242627C8

## Applications

### Tested Applications:

WB, Indirect ELISA

### Species Specificity:

human, mouse

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85340-2-RR (DTX3L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85340-2-PBS in a different storage buffer formulation.