

For Research Use Only

DTX3L Recombinant antibody

Catalog Number: 85340-2-RR



Basic Information

Catalog Number:

85340-2-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2568

GenBank Accession Number:

BC060509

GeneID (NCBI):

151636

UNIPROT ID:

Q8TDB6

Full Name:

deltex 3-like (Drosophila)

Calculated MW:

740 aa, 84 kDa

Observed MW:

84 kDa

Purification Method:

Protein A purification

CloneNo.:

242627C8

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells, HeLa cells, HepG2 cells, MCF-7 cells, U2OS cells, mouse spleen tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85340-2-RR (DTX3L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.