

For Research Use Only

MYH3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 85369-1-PBS



Basic Information

Catalog Number:

85369-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002470

GeneID (NCBI):

4621

UNIPROT ID:

P11055

Full Name:

myosin, heavy chain 3, skeletal muscle, embryonic

Calculated MW:

234 kDa

Observed MW:

220-230 kDa

Purification Method:

Protein A purification

CloneNo.:

242722C1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85369-1-RR (MYH3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85369-1-PBS in a different storage buffer formulation.