

For Research Use Only

MYH3 Recombinant monoclonal antibody

Catalog Number: 85369-1-RR



Basic Information

Catalog Number: 85369-1-RR	GenBank Accession Number: NM_002470	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 4621	CloneNo.: 242722C1
Source: Rabbit	UNIPROT ID: P11055	Recommended Dilutions: WB: 1:500-1:2000
Isotype: IgG	Full Name: myosin, heavy chain 3, skeletal muscle, embryonic	
	Calculated MW: 234 kDa	
	Observed MW: 220-230 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse embryo tissue, mouse eye tissue
Species Specificity: human, mouse	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85369-1-RR (MYH3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.