

For Research Use Only

# CD151 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85385-3-PBS



## Basic Information

<b>Catalog Number:</b> 85385-3-PBS	<b>GenBank Accession Number:</b> NM_004357.5	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 977	<b>CloneNo.:</b> 252763G12
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P48509	
<b>Isotype:</b> IgG	<b>Full Name:</b> CD151 molecule (Raph blood group)	
<b>Immunogen Catalog Number:</b> EG3476	<b>Calculated MW:</b> 28kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
human

## Product Information

85385-3-PBS targets CD151 as part of a matched antibody pair:

MP03215-1: 85385-4-PBS capture and 85385-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean CD151 concentration was determined to be 75.35 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load, 40.21 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load, 128.99 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load, 6.04 ng/mL in Daudi cell extract based on a 1.10 mg/mL extract load, 489.64 ng/mL in A431 cell extract based on a 1.50 mg/mL extract load and 19.76 ng/mL in A549 cell extract based on a 1.20 mg/mL extract



Sandwich ELISA standard curve of MP03215-1, Human CD151 Recombinant Matched Antibody Pair - PBS only. 85385-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3476. 85385-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.