

For Research Use Only

Cas9 Recombinant antibody, PBS Only

Catalog Number: 85389-1-PBS



Basic Information

Catalog Number:

85389-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25281

GenBank Accession Number:

NP_269215

GeneID (NCBI):

Full Name:

Calculated MW:

160 kDa

Purification Method:

Protein A purification

CloneNo.:

242781G12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

n/a

Background Information

Cas9 (CRISPR associated protein 9) is an RNA-guided DNA endonuclease enzyme associated with the CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) adaptive immunity system in *Streptococcus pyogenes*. This antibody detects the sp Cas9.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85389-1-RR (Cas9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85389-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85389-1-RR against Human Cas9 were performed. The affinity constant is 3.18 nM.