

For Research Use Only

AUH Recombinant antibody, PBS Only

Catalog Number: 85421-1-PBS



Basic Information

Catalog Number:

85421-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11265

GenBank Accession Number:

BC020722

GeneID (NCBI):

549

UNIPROT ID:

Q13825

Full Name:

AU RNA binding protein/enoyl-
Coenzyme A hydratase

Calculated MW:

310 aa, 33 kDa

Observed MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

242854F12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85421-1-RR (AUH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85421-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85421-1-RR against Human AUH were performed. The affinity constant is 29.3 nM.