

For Research Use Only

CD158G/KIR2DS5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85436-2-PBS



Basic Information

Catalog Number:

85436-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014513.3

GeneID (NCBI):

3810

UNIPROT ID:

Q14953

Full Name:

killer cell immunoglobulin-like receptor, two domains, short cytoplasmic tail, 5

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

242877F5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

85436-2-PBS targets CD158G/KIR2DS5 as part of a matched antibody pair:

MP01954-1: 85436-1-PBS capture and 85436-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01954-1, CD158G/KIR2DS5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85436-1-PBS. Detection antibody: 85436-2-PBS. Standard: Eg3719. Range: 0.781-100 ng/mL.