

For Research Use Only

PI4K2A Recombinant antibody

Catalog Number: 85448-1-RR



Basic Information

Catalog Number:

85448-1-RR

Size:

100ul , Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7564

GenBank Accession Number:

BC003167

GeneID (NCBI):

55361

UNIPROT ID:

Q9BTU6

Full Name:

phosphatidylinositol 4-kinase type 2

alpha

Calculated MW:

54 kDa

Observed MW:

48-52 kDa

Purification Method:

Protein A purification

CloneNo.:

242992A10

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HeLa cells, HCT 116 cells, HL-60 cells, K-562 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85448-1-RR (PI4K2A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 85448-1-RR against Human PI4K2A were performed. The affinity constant is 2.68 nM.