

For Research Use Only

PI4K2A Recombinant antibody, PBS Only (Capture)

Catalog Number: 85448-5-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 85448-5-PBS | GenBank Accession Number: BC003167 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 55361 | CloneNo.: 242992C6 |
| Source: Rabbit | UNIPROT ID: Q9BTU6 | |
| Isotype: IgG | Full Name: phosphatidylinositol 4-kinase type 2 alpha | |
| Immunogen Catalog Number: AG7564 | Calculated MW: 54 kDa | |

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

85448-5-PBS targets PI4K2A as part of a matched antibody pair:

MP01939-3: 85448-5-PBS capture and 85448-4-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01939-3, Human PI4K2A Recombinant Matched Antibody Pair - PBS only. 85448-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7564. 85448-4-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



The mean PI4K2A concentration was determined to be 158.52 ng/mL in HL-60 cell extract based on a 1.50 mg/mL extract load.