

For Research Use Only

# GLUD2 Recombinant antibody, PBS Only

Catalog Number: 85602-5-PBS



## Basic Information

### Catalog Number:

85602-5-PBS

### Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG5748

### GenBank Accession Number:

BC050732

### GeneID (NCBI):

2747

### UNIPROT ID:

P49448

### Full Name:

glutamate dehydrogenase 2

### Calculated MW:

61 kDa

### Observed MW:

50-60 kDa

### Purification Method:

Protein A purification

### CloneNo.:

243041D11

## Applications

### Tested Applications:

WB, Indirect ELISA

### Species Specificity:

human

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85602-5-RR (GLUD2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85602-5-PBS in a different storage buffer formulation.