

For Research Use Only

BST2 Recombinant antibody

Catalog Number: 85614-4-RR



Basic Information

| | | |
|--|---|---|
| Catalog Number: 85614-4-RR | GenBank Accession Number: NM_198095.2 | Purification Method: Protein A purification |
| Size: 100ul , Concentration: 1000 ug/ml by Nanodrop; | GeneID (NCBI): 69550 | CloneNo.: 242938E6 |
| Source: Rabbit | UNIPROT ID: Q8R2Q8 | Recommended Dilutions: WB: 1:2000-1:10000 |
| Isotype: IgG | Full Name: bone marrow stromal cell antigen 2 | |
| Immunogen Catalog Number: EG2669 | Calculated MW: 19 kDa | |
| | Observed MW: 30-36 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB, ELISA | Positive Controls: |
| Species Specificity: mouse | WB : RAW 264.7 cells, mouse spleen tissue, mouse thymus tissue |

Background Information

BST2, also named as CD317 and Tetherin, belongs to the tetherin family. It may be involved in the sorting of secreted proteins and it is involved in pre-B-cell growth. BST2 is an antiretroviral defense protein, that blocks release of retrovirus from the cell surface. Depleted upon HIV-1 infection by viral VPU protein through 20S proteasome degradation. Depleted upon infection by human Kaposi's sarcoma-associated herpesvirus (KSHV) through ubiquitination and subsequent degradation. BST2 may play a role in B-cell activation in rheumatoid arthritis. It is recently identified interferon-induced cellular proteins that restrict infections by retroviruses and filoviruses and of influenza virus and flaviviruses, respectively. BST2 is a plasma membrane proteins, tetherin inhibits virion particle release from infected cells. BST2 is effective against retroviruses and flavoviruses whilst IFITMs disrupt influenza and flavivirus infection. Observed MW of BST2 is 30-36 kDa (PMID: 19196977; 21237475).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85614-4-RR (Bst2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.