

For Research Use Only

Mouse Cd300lg Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85624-3-PBS



Basic Information

Catalog Number: 85624-3-PBS	GenBank Accession Number: NM_001160712.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 52685	CloneNo.: 242032A9
Source: Rabbit	UNIPROT ID: NP_001154183.1	
Isotype: IgG	Full Name: CD300 antigen like family member G	
Immunogen Catalog Number: EG2731	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

85624-3-PBS targets Cd300lg as part of a matched antibody pair:

MP01979-1: 85624-3-PBS capture and 85624-2-PBS detection (validated in Cytometric bead array)

MP01979-2: 85624-3-PBS capture and 85624-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01979-1, MOUSE Cd300lg Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85450-1-PBS. Detection antibody: 85450-2-PBS. Standard: Eg2731. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01979-2, MOUSE Cd300lg Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85450-1-PBS. Detection antibody: 85450-3-PBS. Standard: Eg2731. Range: 0.625-80 ng/mL