

For Research Use Only

IQGAP1 Recombinant antibody, PBS Only

Catalog Number: 85651-1-PBS



Basic Information

Catalog Number:

85651-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19321

GenBank Accession Number:

BC151834

GeneID (NCBI):

8826

ENSEMBL Gene ID:

ENSG00000140575

UNIPROT ID:

P46940

Full Name:

IQ motif containing GTPase activating protein 1

Calculated MW:

189 kDa

Purification Method:

Protein A purification

CloneNo.:

242206G2

Applications

Tested Applications:

IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunofluorescent analysis of (4% PFA) fixed A431 cells using IQGAP1 antibody (85651-1-RR, Clone: 242206G2) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 85651-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 85651-1-RR (IQGAP1 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 85651-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 85651-1-RR (IQGAP1 antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 85651-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using IQGAP1 antibody (85651-1-RR, Clone: 242206G2) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 85651-1-PBS in a different storage buffer formulation.