

For Research Use Only

PARVB Recombinant antibody, PBS Only

Catalog Number: 85711-1-PBS



Basic Information

Catalog Number:

85711-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5760

GenBank Accession Number:

BC046103

GeneID (NCBI):

29780

UNIPROT ID:

Q9HBI1

Full Name:

parvin, beta

Calculated MW:

42 kDa

Observed MW:

42 kDa, 45 kDa

Purification Method:

Protein A purification

CloneNo.:

242563A9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85711-1-RR (PARVB antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85711-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85711-1-RR against Human PARVB were performed. The affinity constant is below 1 pM.