

For Research Use Only

PARVB Recombinant antibody

Catalog Number: 85711-1-RR



Basic Information

Catalog Number:

85711-1-RR

Size:

100ul , Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5760

GenBank Accession Number:

BC046103

GeneID (NCBI):

29780

UNIPROT ID:

Q9HBI1

Full Name:

parvin, beta

Calculated MW:

42 kDa

Observed MW:

42 kDa, 45 kDa

Purification Method:

Protein A purification

CloneNo.:

242563A9

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, U-87 MG cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85711-1-RR (PARVB antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 85711-1-RR against Human PARVB were performed. The affinity constant is below 1 pM.