

For Research Use Only

SPATA2 Recombinant antibody, PBS Only

Catalog Number: 85712-1-PBS



Basic Information

Catalog Number:

85712-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30773

GenBank Accession Number:

BC009481

GeneID (NCBI):

9825

UNIPROT ID:

Q9UM82

Full Name:

spermatogenesis associated 2

Calculated MW:

58 kDa

Observed MW:

58-60 kDa

Purification Method:

Protein A purification

CloneNo.:

242523H4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85712-1-RR (SPATA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85712-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85712-1-RR against Human SPATA2 were performed. The affinity constant is 0.175 nM.