

For Research Use Only

# HNRNPUL2-Specific Recombinant antibody, PBS Only

Catalog Number: 85739-1-PBS



## Basic Information

### Catalog Number:

85739-1-PBS

### Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### GenBank Accession Number:

NM\_001079559

### GeneID (NCBI):

221092

### UNIPROT ID:

Q1KMD3

### Full Name:

heterogeneous nuclear ribonucleoprotein U-like 2

### Calculated MW:

85 kDa

### Purification Method:

Protein A purification

### CloneNo.:

250029G10

## Applications

### Tested Applications:

IF/ICC, Indirect ELISA

### Species Specificity:

human

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HNRNPUL2-Specific antibody (85739-1-RR, Clone: 250029G10 ) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85739-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using HNRNPUL2-Specific antibody (85739-1-RR, Clone: 250029G10 ) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85739-1-PBS in a different storage buffer formulation.