

For Research Use Only

HNRNPUL2-Specific Recombinant antibody

Catalog Number: 85739-1-RR



Basic Information

Catalog Number:

85739-1-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001079559

GeneID (NCBI):

221092

UNIPROT ID:

Q1KMD3

Full Name:

heterogeneous nuclear ribonucleoprotein U-like 2

Calculated MW:

85 kDa

Purification Method:

Protein A purification

CloneNo.:

250029G10

Recommended Dilutions:

IF/ICC 1:250-1:1000

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells, Jurkat cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HNRNPUL2-Specific antibody (85739-1-RR, Clone: 250029G10) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using HNRNPUL2-Specific antibody (85739-1-RR, Clone: 250029G10) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).