

For Research Use Only

HNRNPUL2-Specific Recombinant antibody, PBS Only

Catalog Number: 85739-2-PBS



Basic Information

Catalog Number:

85739-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001079559

GeneID (NCBI):

221092

UNIPROT ID:

Q1KMD3

Full Name:

heterogeneous nuclear ribonucleoprotein U-like 2

Calculated MW:

85 kDa

Observed MW:

110-120 kDa

Purification Method:

Protein A purification

CloneNo.:

250029E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85739-2-RR (HNRNPUL2-Specific antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85739-2-PBS in a different storage buffer formulation.