

For Research Use Only

RIPK2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85760-5-PBS



Basic Information

Catalog Number: 85760-5-PBS	GenBank Accession Number: BC004553	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8767	CloneNo.: 250069B7
Source: Rabbit	UNIPROT ID: O43353	
Isotype: IgG	Full Name: receptor-interacting serine-threonine kinase 2	
Immunogen Catalog Number: AG0909	Calculated MW: 61 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

85760-5-PBS targets RIPK2 as part of a matched antibody pair:

MP02093-3: 85760-6-PBS capture and 85760-5-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.



Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02093-3, Human RIPK2 Recombinant Matched Antibody Pair - PBS only. 85760-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0909. 85760-5-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL.

The mean RIPK2 concentration was determined to be 76.00 ng/mL in THP-1 cell extract based on a 1.50 mg/mL extract load.