

For Research Use Only

SCN7A-Specific Recombinant antibody, PBS Only

Catalog Number: 85797-1-PBS



Basic Information

Catalog Number:

85797-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002976

GeneID (NCBI):

6332

UNIPROT ID:

Q01118

Full Name:

sodium channel, voltage-gated, type
VII, alpha

Calculated MW:

193 kDa

Observed MW:

180 kDa

Purification Method:

Protein A purification

CloneNo.:

250004D8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SCN7A, also named as SCN6A, belongs to the sodium channel family. SCN7A mediates the voltage-dependent sodium ion permeability of excitable membranes. Assuming opened or closed conformations in response to the voltage difference across the membrane, SCN7A forms a sodium-selective channel through which Na⁺ ions may pass in accordance with their electrochemical gradient. This antibody is specific to SCN7A.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 85797-1-RR (SCN7A-Specific antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85797-1-PBS in a different storage buffer formulation.



mouse heart tissue were subjected to SDS PAGE followed by western blot with 85797-1-RR (SCN7A-Specific antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85797-1-PBS in a different storage buffer formulation.