

For Research Use Only

# MEF2D Recombinant antibody, PBS Only (Detector)

Catalog Number: 85801-2-PBS



## Basic Information

<b>Catalog Number:</b> 85801-2-PBS	<b>GenBank Accession Number:</b> BC040949	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 4209	<b>CloneNo.:</b> 250068G12
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q14814	
<b>Isotype:</b> IgG	<b>Full Name:</b> myocyte enhancer factor 2D	
<b>Immunogen Catalog Number:</b> AG5700	<b>Calculated MW:</b> 56 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

85801-2-PBS targets MEF2D as part of a matched antibody pair:

MP02136-1: 85801-1-PBS capture and 85801-2-PBS detection (validated in Cytometric bead array)

MP02136-2: 85801-3-PBS capture and 85801-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MPO2136-2, MEF2D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85801-3-PBS. Detection antibody: 85801-2-PBS. Standard: Ag5700. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MPO2136-1, MEF2D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85801-1-PBS. Detection antibody: 85801-2-PBS. Standard: Ag5700. Range: 0.781-100 ng/mL.