

For Research Use Only

MOV10 Recombinant antibody, PBS Only

Catalog Number: 85818-2-PBS



Basic Information

Catalog Number:

85818-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0384

GenBank Accession Number:

BC002548

GeneID (NCBI):

4343

UNIPROT ID:

Q9HCE1

Full Name:

Mov10, Moloney leukemia virus 10, homolog (mouse)

Calculated MW:

114 kDa

Purification Method:

Protein A purification

CloneNo.:

250130F7

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

