

For Research Use Only

C6orf15 Recombinant antibody

Catalog Number: 85830-1-RR



Basic Information

Catalog Number: 85830-1-RR	GenBank Accession Number: BC153157	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 29113	CloneNo.: 250126F4
Source: Rabbit	UNIPROT ID: Q6UXA7	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG	Full Name: chromosome 6 open reading frame 15	
Immunogen Catalog Number: AG20276	Calculated MW: 325 aa, 34 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HaCaT cells, Caco-2 cells, PC-3 cells
Species Specificity: human	

Background Information

C6orf15, also named as Protein STG, is a 325 amino acid protein, which as a secreted protein is expressed in skin and tonsils. The calculated molecular weight of C6orf15 is 34 kDa, but we detected a 50-55 kDa protein by western blot. C6orf15 may contains some modified residues.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85830-1-RR (C6orf15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.