

For Research Use Only

Rat Fas Ligand/TNFSF6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85836-3-PBS



Basic Information

Catalog Number: 85836-3-PBS	GenBank Accession Number: NM_012908.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 25385	CloneNo.: 242892H8
Source: Rabbit	UNIPROT ID: P36940	
Isotype: IgG	Full Name: Fas ligand (TNF superfamily, member 6)	
Immunogen Catalog Number: EG3218	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Product Information

85836-3-PBS targets Fas Ligand/TNFSF6 as part of a matched antibody pair:

MP02135-1: 85836-4-PBS capture and 85836-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2135-1, RAT Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85836-4-PBS. Detection antibody: 85836-3-PBS. Standard: Eg3218. Range: 1.563-200 ng/mL.