

For Research Use Only

Mouse Plat Recombinant antibody, PBS Only

Catalog Number:85872-1-PBS



Basic Information

Catalog Number: 85872-1-PBS	GenBank Accession Number: NM_008872.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 18791	CloneNo.: 243116G3
Source: Rabbit	UNIPROT ID: P11214	
Isotype: IgG	Full Name: plasminogen activator, tissue	
Immunogen Catalog Number: EG3025	Calculated MW: 63 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
mouse

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data