

For Research Use Only

LGALS3BP Recombinant antibody, PBS Only (Detector)

Catalog Number: 85883-1-PBS



Basic Information

Catalog Number: 85883-1-PBS	GenBank Accession Number: NM_005567	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3959	CloneNo.: 250138D9
Source: Rabbit	UNIPROT ID: Q08380	
Isotype: IgG	Full Name: lectin, galactoside-binding, soluble, 3 binding protein	
Immunogen Catalog Number: EG3302	Calculated MW: 585 aa, 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85883-1-PBS targets LGALS3BP as part of a matched antibody pair:

MP02172-1: 85883-2-PBS capture and 85883-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2172-1, LGALS3BP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85883-2-PBS. Detection antibody: 85883-1-PBS. Standard: Eg3302. Range: 0.781-25 ng/mL