

For Research Use Only

# TAF1B Recombinant antibody, PBS Only (Capture)

Catalog Number: 85885-3-PBS



## Basic Information

### Catalog Number:

85885-3-PBS

### Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG38211

### GenBank Accession Number:

BC018137

### GeneID (NCBI):

9014

### UNIPROT ID:

Q53T94

### Full Name:

TATA box binding protein (TBP)-associated factor, RNA polymerase I, B, 63kDa

### Calculated MW:

588 aa, 69 kDa

### Purification Method:

Protein A purification

### CloneNo.:

250162B8

## Applications

### Tested Applications:

Cytometric bead array, Indirect ELISA

### Species Specificity:

human

## Product Information

85885-3-PBS targets TAF1B as part of a matched antibody pair:

MP02183-2: 85885-3-PBS capture and 85885-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

### Storage:

Store at -80°C.

### Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02183-2, TAF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85885-3-PBS. Detection antibody: 85885-2-PBS. Standard: Ag38211. Range: 0.781-100 ng/mL.