

For Research Use Only

SNX15 Recombinant antibody, PBS Only

Catalog Number: 85909-3-PBS



Basic Information

Catalog Number:

85909-3-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9006

GenBank Accession Number:

BC009897

GeneID (NCBI):

29907

UNIPROT ID:

Q9NRS6

Full Name:

sorting nexin 15

Calculated MW:

342 aa, 38 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

250152G10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

SNX15 belongs to the sorting nexin family and is involved in several stages of intracellular trafficking. Overexpression of SNX15 results in a decrease in the processing of insulin and hepatocyte growth factor receptors to their mature subunits (PMID: 11085978). SNX15 is widely expressed, with the highest expression in skeletal muscle, heart, brain, kidney, spleen, thymus, and small intestine.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85909-3-RR (SNX15 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85909-3-PBS in a different storage buffer formulation.