

For Research Use Only

DGUOK Recombinant antibody, PBS Only

Catalog Number: 85932-2-PBS



Basic Information

Catalog Number:

85932-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9361

GenBank Accession Number:

BC015757

GeneID (NCBI):

1716

UNIPROT ID:

Q16854

Full Name:

deoxyguanosine kinase

Calculated MW:

277 aa, 32 kDa

Observed MW:

28-32 kDa

Purification Method:

Protein A purification

CloneNo.:

250161D11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85932-2-RR (DGUOK antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85932-2-PBS in a different storage buffer formulation.