

For Research Use Only

MGAT2 Recombinant antibody, PBS Only

Catalog Number: 85933-1-PBS



Basic Information

Catalog Number:

85933-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8423

GenBank Accession Number:

BC006390

GeneID (NCBI):

4247

UNIPROT ID:

Q10469

Full Name:

mannosyl (alpha-1,6-)-glycoprotein beta-1,2-N-acetylglucosaminyltransferase

Calculated MW:

447 aa, 52 kDa

Observed MW:

55-60 kDa

Purification Method:

Protein A purification

CloneNo.:

250147B1

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85933-1-RR (MGAT2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85933-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Raji cells using MGAT2 antibody (85933-1-RR, Clone: 250147B1) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 85933-1-PBS in a different storage buffer formulation.