

For Research Use Only

SSB Recombinant antibody, PBS Only

Catalog Number: 85944-2-PBS



Basic Information

Catalog Number:

85944-2-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2294

GenBank Accession Number:

BC001289

GeneID (NCBI):

6741

UNIPROT ID:

P05455

Full Name:

Sjogren syndrome antigen B (autoantigen La)

Calculated MW:

408 aa, 47 kDa

Observed MW:

46-50 kDa

Purification Method:

Protein A purification

CloneNo.:

250465G10

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using SSB antibody (85944-2-RR, Clone: 250465G10) at dilution of 1:1000 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85944-2-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 85944-2-RR (SSB antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85944-2-PBS in a different storage buffer formulation.