

For Research Use Only

# IL-23 Recombinant antibody

Catalog Number: 85947-1-RR



## Basic Information

Catalog Number:

85947-1-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1254

GenBank Accession Number:

GeneID (NCBI):

51561

UNIPROT ID:

Q9NPF7

Full Name:

interleukin 23, alpha subunit p19

Observed MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

250414B6

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein,

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 85947-1-RR (IL-23 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.