

For Research Use Only

IL-23 Recombinant antibody, PBS Only

Catalog Number: 85947-2-PBS



Basic Information

Catalog Number:

85947-2-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1254

GenBank Accession Number:

GeneID (NCBI):

51561

UNIPROT ID:

Q9NPF7

Full Name:

interleukin 23, alpha subunit p19

Purification Method:

Protein A purification

CloneNo.:

250414G4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85947-2-RR (IL-23 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85947-2-PBS in a different storage buffer formulation.



Recombinant protein were subjected to SDS PAGE followed by western blot with 85947-2-RR (IL-23 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85947-2-PBS in a different storage buffer formulation.