

For Research Use Only

IL-23 p19 Recombinant monoclonal antibody

Catalog Number: 85947-3-RR



Basic Information

Catalog Number: 85947-3-RR	GenBank Accession Number: GeneID (NCBI): 51561	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	UNIPROT ID: Q9NPF7	CloneNo.: 252904G4
Source: Rabbit	Full Name: interleukin 23, alpha subunit p19	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG	Observed MW: 21 kDa	
Immunogen Catalog Number: HZ-1254		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : PMA, LPS and Brefeldin A treated THP-1 cells,
Species Specificity: human	

Background Information

Interleukin 23 (IL-23) is a member of the IL12 cytokine family and composed of two subunits, IL12p40 and IL23p19. It is produced by antigen presenting cells and has been shown to promote the production and survival of a distinct lineage of T-cells called TH17 cells. A functional receptor for IL-23 (the IL-23 receptor) has been identified and is composed of IL-12Rβ1 and IL-23R. IL-23 is expressed chiefly by the macrophages and DCs. The IL-23R is found on memory T cells, NKT cells, macrophages, DCs, and naive T cells upon activation by TGF-β and IL-6. The main biological effects of IL-23 identified initially consist of stimulation of antigen presentation by DCs, T cell differentiation to Th17 cells, and production of interferon-γ (IFN-γ). IL-23 acts also as an end-stage effector cytokine through direct action on macrophages.

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PMA, LPS and Brefeldin A treated THP-1 cells were subjected to SDS PAGE followed by western blot with 85947-3-RR (IL-23 p19 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.