

For Research Use Only

SLC22A15 Recombinant antibody, PBS Only

Catalog Number: 85960-1-PBS



Basic Information

Catalog Number:

85960-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14547

GenBank Accession Number:

BC026358

GeneID (NCBI):

55356

UNIPROT ID:

Q8IZD6

Full Name:

solute carrier family 22, member 15

Calculated MW:

535aa,59 kDa; 547aa,61 kDa

Observed MW:

61-66 kDa

Purification Method:

Protein A purification

CloneNo.:

250347G3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85960-1-RR (SLC22A15 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85960-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85960-1-RR against Human SLC22A15 were performed. The affinity constant is 5.00 nM.