

For Research Use Only

PAPP A2 Recombinant antibody

Catalog Number: 85976-3-RR



Basic Information

Catalog Number:

85976-3-RR

Size:

100ul, Concentration: 1000 µg/ml by 60676

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2812

GenBank Accession Number:

NM_020318.2

GeneID (NCBI):

60676

UNIPROT ID:

Q9BXP8-1

Full Name:

pappalysin 2

Calculated MW:

199 kDa

Observed MW:

245 kDa

Purification Method:

Protein A purification

CloneNo.:

250602C12

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue, fetal human brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 85976-3-RR (PAPPA2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.