

For Research Use Only

ALDH1L1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85980-4-PBS



Basic Information

Catalog Number: 85980-4-PBS	GenBank Accession Number: BC027241	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10840	CloneNo.: 250565H6
Source: Rabbit	UNIPROT ID: O75891	
Isotype: IgG	Full Name: aldehyde dehydrogenase 1 family, member L1	
Immunogen Catalog Number: AG11285	Calculated MW: 99 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85980-4-PBS targets ALDH1L1 as part of a matched antibody pair:

MP02776-1: 85980-5-PBS capture and 85980-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ALDH1L1 concentration was determined to be 4047.1 ng/mL in Mouse liver tissue extract based on a 8.4 mg/mL extract load, 103.8 ng/mL in Mouse brain tissue extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02776-1, ALDH1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85980-5-PBS. Detection antibody: 85980-4-PBS. Standard: Ag11285. Range: 0.781-100 ng/mL