

For Research Use Only

BRD2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85990-3-PBS



Basic Information

Catalog Number: 85990-3-PBS	GenBank Accession Number: BC063840	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6046	CloneNo.: 250453A12
Source: Rabbit	UNIPROT ID: P25440	
Isotype: IgG	Full Name: bromodomain containing 2	
Immunogen Catalog Number: AG17589	Calculated MW: 836 aa, 92 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85990-3-PBS targets BRD2 as part of a matched antibody pair:

MP02788-1: 85990-4-PBS capture and 85990-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BRD2 concentration was determined to be 23.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 47.6 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 21.7 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02788-1, BRD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85990-4-PBS. Detection antibody: 85990-3-PBS. Standard: Ag17589. Range: 1.25-40 ng/mL