

For Research Use Only

Mouse ANGPTL4 Recombinant antibody, PBS Only

Catalog Number: 86006-4-PBS



Basic Information

Catalog Number: 86006-4-PBS	GenBank Accession Number: NM_020581.2	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 57875	CloneNo.: 250335F5
Source: Rabbit	UNIPROT ID: Q9Z1P8	
Isotype: IgG	Full Name: angiopoietin-like 4	
Immunogen Catalog Number: EG2770	Calculated MW: 46 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
mouse, rat

Background Information

ANGPTL4 was first identified in adipose tissue and later found to be widely expressed in several tissues including the liver, skeletal muscle, small intestine, placenta, heart, and kidney. ANGPTL4 expression was first identified to be regulated by peroxisome proliferator-activator receptors (PPARs), PPAR α and PPAR γ .

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86006-4-RR (Angptl4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86006-4-PBS in a different storage buffer formulation.