

For Research Use Only

COPZ1 Recombinant antibody, PBS Only

Catalog Number: 86014-1-PBS



Basic Information

Catalog Number:

86014-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14234

GenBank Accession Number:

BC002849

GeneID (NCBI):

22818

UNIPROT ID:

P61923

Full Name:

coatamer protein complex, subunit zeta 1

Calculated MW:

177 aa, 20 kDa

Observed MW:

15-20 kDa

Purification Method:

Protein A purification

CloneNo.:

250388F3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, monkey

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86014-1-RR (COPZ1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86014-1-PBS in a different storage buffer formulation.