

For Research Use Only

COPZ1 Recombinant antibody

Catalog Number: 86014-1-RR



Basic Information

Catalog Number:

86014-1-RR

GenBank Accession Number:

BC002849

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

GeneID (NCBI):

22818

CloneNo.:

250388F3

Source:

Rabbit

UNIPROT ID:

P61923

Recommended Dilutions:

WB: 1:5000-1:50000

Isotype:

IgG

Full Name:

coatomer protein complex, subunit zeta 1

Immunogen Catalog Number:

AG14234

Calculated MW:

177 aa, 20 kDa

Observed MW:

15-20 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : HeLa cells, RAW 264.7 cells, COS-7 cells

Species Specificity:

human, mouse, monkey

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86014-1-RR (COPZ1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.