

For Research Use Only

DIS3L2 Recombinant antibody, PBS Only

Catalog Number: 86034-1-PBS



Basic Information

Catalog Number:

86034-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30634

GenBank Accession Number:

BC036113

GeneID (NCBI):

129563

UNIPROT ID:

Q8IYB7

Full Name:

DIS3 mitotic control homolog (S. cerevisiae)-like 2

Observed MW:

99 kDa

Purification Method:

Protein A purification

CloneNo.:

250578A2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86034-1-RR (DIS3L2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86034-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86034-1-RR against Human DIS3L2 were performed. The affinity constant is 0.463 nM.