

For Research Use Only

# DIS3L2 Recombinant antibody

Catalog Number: 86034-1-RR



## Basic Information

Catalog Number:

86034-1-RR

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30634

GenBank Accession Number:

BC036113

GeneID (NCBI):

129563

UNIPROT ID:

Q8IYB7

Full Name:

DIS3 mitotic control homolog (S. cerevisiae)-like 2

Observed MW:

99 kDa

Purification Method:

Protein A purification

CloneNo.:

250578A2

Recommended Dilutions:

WB: 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, K-562 cells, A2780 cells

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86034-1-RR (DIS3L2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86034-1-RR against Human DIS3L2 were performed. The affinity constant is 0.463 nM.