

For Research Use Only

FOXJ1 Recombinant antibody

Catalog Number: 86070-1-RR



Basic Information

Catalog Number:

86070-1-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28294

GenBank Accession Number:

BC046460

GeneID (NCBI):

2302

UNIPROT ID:

Q92949

Full Name:

forkhead box J1

Calculated MW:

45 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

250724E9

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, Caco-2 cells, T-47D cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86070-1-RR (FOXJ1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.