

For Research Use Only

MAG/Siglec-4a Recombinant antibody

Catalog Number: 86073-1-RR



Basic Information

Catalog Number:

86073-1-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3251

GenBank Accession Number:

NM_017190.4

GeneID (NCBI):

29409

UNIPROT ID:

P07722-1

Full Name:

myelin-associated glycoprotein

Calculated MW:

69 kDa

Observed MW:

70-100 kDa

Purification Method:

Protein A purification

CloneNo.:

250740D11

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : rat brain tissue, mouse brain tissue, rat spinal cord tissue, mouse spinal cord tissue, rat cerebellum tissue, mouse cerebellum tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86073-1-RR (MAG/Siglec-4a antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.