

For Research Use Only

Mouse Ccl8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86081-2-PBS



Basic Information

Catalog Number: 86081-2-PBS	GenBank Accession Number: NM_021443.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 20307	CloneNo.: 250174H8
Source: Rabbit	UNIPROT ID: Q9Z121	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 8	
Immunogen Catalog Number: EG2771	Calculated MW: 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Product Information

86081-2-PBS targets Ccl8 as part of a matched antibody pair.

MP02245-1: 86081-2-PBS capture and 86081-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Ccl8 concentration was determined to be 1317.1 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02245-1, MOUSE Ccl8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86081-2-PBS. Detection antibody: 86081-1-PBS. Standard: Eg2771. Range: 0.156-20 ng/mL.