

For Research Use Only

PRPF6 Recombinant antibody, PBS Only

Catalog Number: 86123-1-PBS



Basic Information

Catalog Number:

86123-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13817

GenBank Accession Number:

BC001666

GeneID (NCBI):

24148

UNIPROT ID:

O94906

Full Name:

PRP6 pre-mRNA processing factor 6 homolog (S. cerevisiae)

Calculated MW:

941 aa, 107 kDa

Observed MW:

100-108 kDa

Purification Method:

Protein A purification

CloneNo.:

250706F2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

PRPF6 is involved in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex, it can enhance dihydrotestosterone-induced transactivation activity of AR, as well as dexamethasone-induced transactivation activity of NR3C1, but does not affect estrogen-induced transactivation.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86123-1-RR (PRPF6 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86123-1-PBS in a different storage buffer formulation.